

Leakage Tester GB-M2



Introduction

The GB-M2 sealing tester uses compressed air as the power source and high-precision electronic pressure sensor as the measurement method. Through microcomputer control, it can effectively compare and evaluate the sealing process and sealing performance of flexible packages, and provide scientific information for determining relevant technical requirements. basis. It is suitable for sealing performance test of flexible packages in food, pharmaceutical, daily chemical and other industries.

Test principle

By evacuating the vacuum chamber, the sample immersed in water will generate a pressure difference between the inside and outside, and observe the gas escape or water infiltration in the sample, so as to judge the sealing performance of the sample.

Standard

ASTM D3078

Specification

Item	Technical Parameters
Vacuum degree	—(0~90)KPa
Vacuum precision	±1%
Vacuum keeping	1 s ~ 99 h59 m59 s
Vacuum chamber available size	Standard (Φ270*270) mm



	Standard (Φ300*300) mm,
Vacuum tank size	Other sizes optional or customized: middle Φ450*350, large Φ600*500
Intake pressure	(0.6~0.8)MPa
Power supply	220 V, 50 Hz
Instrument size	(317×305×170) mm
Power	220V, 50Hz

Features

◆ High-quality accessories, quality assurance

The built-in high-precision pressure sensor has high sensitivity and can monitor the intake pressure in real time.

◆ Fully automatic operating system

Set the vacuum degree and holding time as needed, the instrument is fully automatic for testing, automatic vacuuming, holding the set time, automatic stop and inflation return.

The new air circuit system uses compressed air and is pneumatically driven, which can realize quick inflation and return to the position, and the internal and external air pressures are consistent.

◆ Excellent appearance and convenient operation

The host is embedded with a 7-inch high-resolution color LCD screen, with clear views, sensitive touch controls, and easy operation.

The instrument has its own operating system, which can be tested independently from the computer, and the data can be processed automatically.



Application

	Paper-plastic composite packaging	Sealability test for paper-plastic composite bags and boxes.
	Plastic packaging	Integral tightness test for plastic bottles, tubes, cans and boxes.

Factory configuration

No.	Name	Model NO./Specification	Qty	Unit	Remark
1	Leakage tester	GB-M2	1	pc	
2	Vacuum tank	/	1	pc	
3	Power line	10A 1.8 mts	1	pc	
4	PU pipe	Φ8*Φ5	5	mts	
5	Reducer Straight Fitting	PG8-6	1	pc	

Note: Guangzhou Biaoji has always been committed to the innovation and improvement of product performance and function. For this reason, product technical specifications and appearance will also be changed accordingly. The above situation will not be notified. The company reserves the right of modification and final interpretation.





Service Commitment

- ◆ Free packaging inspection technical support.
- ◆ One-year warranty for instrument, free software upgrade service, and lifetime free technical support.
- ◆ Free to the factory to train the operation of the instrument, qualified people issued graduation certificates, free accommodation and food.
- ◆ Instrument “trade-in” policy and provision of backup instruments in case of failure.
- ◆ There is a testing center (CNAS L8185) to provide incoming sample testing and sample data comparison services.
- ◆ Set up a Reference Standards Development Center to provide calibration services.

Company Introduction

Guangzhou Biaoji Packaging Equipment Co., Ltd. was established in 2002 and located in Guangzhou Economic Development Zone. It is a company integrating R&D, production and sales of packaging testing instruments, packaging equipment, modified atmosphere preservation equipment and other products, as well as providing third-party testing and standard material proficiency testing. It is committed to providing comprehensive, professional and high-quality products and technical services for packaging, food, medicine, testing and other industries.

Write industry standards with quality, and become the leader and leader in the field of global packaging testing. The 80,000-mile journey is sailing, and Guangzhou Biaoji will escort you!



Address: No. 1, Minghua 3rd Street, Jinxiu Road,
Guangzhou Economic and Technological Development Zone

Tel: 020-8615 3789, 8615 3790

After-sales hotline: 4007886855

Website: <https://www.gbpitester.com/>

Email: info@gbtest.cn



Website QR code



WeChatQR code

