

Heat Shrink Tester GBK-1





Introduction

It is used to determination of the degree of unrestrained linear thermal shrinkage and dimensional stability at different temperatures. Widely used in industries of plastic material, plastic products, petroleum, chemicals, research institute, inspection agency and so on.

Standard

ASTM D2732, GB/T 13519

Specification

Item	Technical parameters	
Sample size	≤140×140mm	
Temperature range	Room temperature~250°C	
Temperature control accuracy	±0.1°C	
Instrument size	340×400×220mm	
Weight	13.5kg	
Power	1500W	
Power supply	AC 220V, 50Hz	

Features

- ◆ High-quality aluminum alloy liquid tank and mask, good heat conduction and heat preservation effect.
- ◆ LCD display, display temperature and timing, convenient parameter setting and simple operation.







Product List

1. Spa	are parts				
No.]	Item	Remark	Quantity	
1	Power line		3*0.75 10A 1.8mts	1pc	
2	Cover		FS-1	1 pc	
3	Heat shrinkable fixtures			1 set	
4	Torx handle			2 pcs	
5	Inner hexagon			1 pc	
6	Phillips screwdriver		6-10mm	1 pc	
7	100x100 sample plate			1 pc	
2. Op	tional part				
No.	Item	Remark	Quantity		
1	Metrology certificate			Optional	
3. Use	er provided		,		
1	Power 500w 220V 10A, 3-hole socket, Ground wire				
2	Silicone oil or glycerin				

Note: GBPI has always been committed to the innovation and improvement of product performance and function. For this reason, product technical specifications and appearance will also be changed accordingly.

The above situation will not be notified. GBPI reserves the right of modification and final interpretation.









Service Commitment

- Free packaging inspection technical support.
- One-year warranty for instrument, free software upgrade service, and lifetime free technical support.
- Free to the factory to train the operation of the instrument, qualified people issued graduation certificates, free accommodation and food.
- ◆ Instrument "trade-in" policy and provision of backup instruments in case of failure.
- There is a testing center (CNAS L8185) to provide incoming sample testing and sample data comparison services.
- Set up a Reference Standards Development
 Center to provide calibration services.

Company Introduction

Guangzhou Biaoji Packaging Equipment Co.,
Ltd. was established in 2002 and located in
Guangzhou Economic Development Zone. It is a
company integrating R&D, production and sales
of packaging testing instruments, packaging
equipment, modified atmosphere preservation
equipment and other products, as well as
providing third-party testing and standard
material proficiency testing. It is committed to
providing comprehensive, professional and
high-quality products and technical services for
packaging, food, medicine, testing and other
industries.

Write industry standards with quality, and become the leader and leader in the field of global packaging testing. The 80,000-mile journey is sailing, and GBPI will escort you!







Address: No. 3, Linjiang Road, Huangpu District, Guangzhou, Guangdong Province, China

Tel.: 86-20-86153790

After Service: 86-136 6047 2664 (Moblie/WhatsApp)

Website: www.gbpitester.com

Email: info@gbtest.cn





Website QR code

WeChat QR code



