

# **Electronic Tensile Tester GBL-H**





#### **Product introduction**

The electronic tensile testing machine GBL-H is a mechanical performance testing machine developed and manufactured by the Guangzhou Biaoji R&D team based on GB, JIS, ASTM and other standard requirements and market demand. The instrument uses high-precision sensors, with force measurement accuracy within ±1% FS; professional software support, fully automatic operation, equipped with an LCD screen, real-time display of test curve data and queryable test reports, making testing convenient and fast; the control measurement unit can Choose built-in or external type, with complete safety protection functions such as limit protection, overload protection, and emergency stop.

It is suitable for stretching, peeling, heat sealing, tearing, etc. of various plastic films, paper-aluminum-plastic composite films, various soft elastic materials, rubber sheets, adhesive tapes, adhesives, paper, non-woven fabrics and other packaging materials. Puncture and other mechanical property tests comply with the testing technology requirements of the national technical supervision department and are widely used in quality inspection, drug inspection, scientific research, packaging, film, food, pharmaceuticals, daily chemicals and other industries.


# **Test principle**

Fix the processed sample between the two chucks of the fixture, start the instrument to make the two chucks move relative to each other, and the force sensor located on the movable chuck collects the force value changes during the test, and the displacement sensor built into the machine The displacement changes are collected to calculate the tensile, peeling, heat sealing, tearing, puncture and other performance indicators of the sample.









Tensile test principle diagram

#### **Standards**

ASTM D828, ASTM E4, ASTM D882, ASTM D1938, ASTM D3330, ASTM F904, JIS P8113, QB/T 2358, GB 8808, GB 13022, GB/T 1040, GB 4850, vGB/T 7753, GB/T 7754, GB/T 453, GB/T 17200, GB/T 16578.1, QB/T 1130, GB/T 2791, GB/T 2790, GB/T 2792, GB/T 7122, GB/T 10004, GB/T 17590, JJG 139, GB/T 6344, GB 10808, YBB 00112003, YBB 00102003, YBB 00132002, YBB 00202004

# **Specification**

Item	Technical Parameters
Force measuring range	0~300 N (0~500 N optional ), 1250% elongation
Force measurement	Indication value is within ±1%FS (Level 1)
accuracy	
Test speed	1-500 mm/min
Displacement accuracy	Indication value is within ±1%
Effective test width	30 mm (50 mm can be customized)
Effective test itinerary	750 mm (1000mm customized)
Size	L×W×H: 475 mm×460 mm×1430 mm
Weight	76.5 kg







Power supply	AC 220 V, 50 Hz
Power	300 W

#### **Features**

### Beautiful and convenient body design

The table-type column structure is beautiful and exquisite, and the seat operation is easy and natural.

#### Flexible and efficient testing methods

It can perform tensile, peeling, heat sealing, tearing, puncture and other tests; multi-parameter variable settings of width, speed, thickness and clamping distance; high sampling rate (50 times/second), the test data is accurate and reliable.

## Safe and high-end mechanical structure

Three-level safety protection; equipped with emergency stop switch, safety is guaranteed; imported core components are used, with high precision and long service life; LCD screen displays test data and working status.

#### **Intelligent operating system**

It can realize user management, authority management and audit tracking, and meets the requirements of GMP Appendix "Computerized System".

The test process parameters can be set flexibly to meet the requirements of non-standard tests and scientific research tests.

It can be tested independently without the computer, the data is automatically processed, and the built-in micro printer can print data in real time. The instrument is equipped with a computer interface and can also be connected to a computer for online testing.

#### **Laboratory intelligent IoT platform**

Supports access to the laboratory intelligent IoT platform, which can realize functions such as managing experimental data, remote diagnosis and troubleshooting.







# Application

Plastic film	It is suitable for tensile, peeling, heat sealing, tearing, puncturing, compression, bending and other performance tests of plastic films, composite materials, flexible packaging materials, tapes and other products.
Composite film bag	It is suitable for tensile, peeling, heat sealing, tearing, puncture, compression, bending and other performance tests of composite films, co-extruded films, composite bags, retort bags, aluminum-plastic composite bags and other products.
Paper	It is suitable for tensile strength, peeling, deformation, tearing, heat sealing, bonding and other performance tests of paper, kraft paper, facial tissue, toilet paper and other products.
Medical Patches	It is suitable for tensile strength, peeling, deformation, tearing, heat sealing, bonding and other performance tests of medical cold compress, plaster and other products.
Pressure Sensitive Tapes	It is suitable for tensile strength, peeling, deformation, tearing, heat sealing, bonding and other performance tests of cellophane tape, electrical insulating tape, medical zinc oxide adhesive tape and other products.







# Configuration

Standard	Power cord, data cable, stretching clamp, puncture clamp, software disk,
configuration	sampling plate, inner hexagon, right-angle tearing knife, Phillips screwdriver
Optional	Computer (Windows system), pneumatic clamp, measurement certificate
parts	
Prepared by	Power supply, 3-hole socket, ground wire, computer (Windows system)
user	

Note: Guangzhou Biaoji has always been committed to the innovation and improvement of product performance and function. For this reason, product technical specifications and appearance will also be changed accordingly. The above situation will not be notified. The company reserves the right of modification and final interpretation.









# **Service Commitment**

- Free packaging inspection technical support.
- One-year warranty for instrument, free software upgrade service, and lifetime free technical support.
- ◆ Free to the factory to train the operation of the instrument, qualified people issued graduation certificates, free accommodation and food.
- Instrument "trade-in" policy and provision of backup instruments in case of failure.
- There is a testing center (CNAS L8185) to provide incoming sample testing and sample data comparison services.
- Set up a Reference Standards Development Center to provide calibration services.

# **Company Introduction**

Guangzhou Biaoji Packaging Equipment Co.,
Ltd. was established in 2002 and located in
Guangzhou Economic Development Zone. It is a
company integrating R&D, production and
sales of packaging testing instruments,
packaging equipment, modified atmosphere
preservation equipment and other products, as
well as providing third-party testing and standard
material proficiency testing. It is committed to
providing comprehensive, professional and
high-quality products and technical services for
packaging, food, medicine, testing and
other industries.

Write industry standards with quality, and become the leader and leader in the field of global packaging testing. The 80, 000-mile journey is sailing, and Guangzhou Biaoji will escort you!







Address: No. 1, Minghua 3rd Street, Jinxiu Road,

Guangzhou Economic and Technological Development Zone

Tel: 020-8615 3789, 8615 3790

After-sales hotline: 4007886855

Website:https://www.gbpitester.com/

Email: info@gbtest.cn





Website QR code

WeChatQR code



