

Sample Cutter for Tensile Tester LC-2



Product Description

The LC-2 Sample Cutter is used for cutting samples of various plastic films, paper and other materials prior to stretching, peeling and heat sealing tests. It is easy and efficient to operate and can accurately cut samples of standard required sizes.

Features

- ◆ Meets the sample making requirements for Type 2 specimens of standard GB/T 1040.3.
- ◆ Using imported tungsten steel finely ground special blades, sharp and durable, with hardness up to 92 HRA.
- ◆ Accurate sample size, no damage, smooth and burr-free edges.
- ◆ Spring tensioned for easy sample fixing and smooth sample cutting.
- ◆ The ability to cut up to 10 samples in a single pass allows for efficient production of large batches of samples.
- ◆ Safety protection, the outer layer of the blade is covered with an aricotta stopper so that the blade does not come into direct contact with it.
- ◆ Different sample widths can be ordered to suit the requirements of each customer.



Specifications

| Item | Specifications |
|-------------------|-------------------------------|
| Cut size | 15 mm×230 mm (Optional width) |
| Number of samples | 10 pcs (Single Layer) |
| Sample thickness | ≤ 200 μm |
| Product size | 250 mm×350 mm×248 mm |
| Weight | 11 kg |

Note: GBPI is always committed to product innovation and improved performance, so accordingly product technical specifications are subject to change without notice. GBPI reserves the right to amend and the final power of interpretation.





Service Commitment

- ◆ Free packaging inspection technical support.
- ◆ One-year warranty for instrument, free software upgrade service, and lifetime free technical support.
- ◆ Free to the factory to train the operation of the instrument, qualified people issued graduation certificates, free accommodation and food.
- ◆ Instrument “trade-in” policy and provision of backup instruments in case of failure.
- ◆ There is a testing center (CNAS L8185) to provide incoming sample testing and sample data comparison services.
- ◆ Set up a Reference Standards Development Center to provide calibration services.

Company Introduction

Guangzhou Biaoji Packaging Equipment Co., Ltd. was established in 2002 and located in Guangzhou Economic Development Zone. It is a company integrating R&D, production and sales of packaging testing instruments, packaging equipment, modified atmosphere preservation equipment and other products, as well as providing third-party testing and standard material proficiency testing. It is committed to providing comprehensive, professional and high-quality products and technical services for packaging, food, medicine, testing and other industries.

Write industry standards with quality, and become the leader and leader in the field of global packaging testing. The 80,000-mile journey is sailing, and GBPI will escort you!



Address: No. 3, Linjiang Road, Huangpu District, Guangzhou, Guangdong Province, China

Tel.: 86-20-86153790

After Service: 86-136 6047 2664 (Mobile/WhatsApp)

Website: www.gbpitester.com

Email: info@gbtest.cn



Website QR code



WeChat QR code

